

R1 shock-proof

Exact precision that lasts! Made in Germany.



Handheld, portable, aneroid sphygmomanometer with shock-proof technology. Approved by the British and Irish Hypertension Society (BIHS). Available with two face scale color options - black or white - with a comprehensive selection of cuff sizes.

- Shock-proof technology to height of 120 cm [47 inches]
- Air duct stability, to ensure precise needle movement
- Fast, ergonomic inflation due to optimized bulb design with integrated spoon
- Non-resistant opening of air release valve ensures precise and comfortable operation
- Scale: Ø 52 mm white with blue pointer or black with red fluorescent pointer
- Metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomic design
- Maximum error tolerance of ± 3 mmHg
- Extensive selection of various cuff sizes
- Supplied in a vinyl zip case

1250-100	R1 shock-proof ® black manometer, 1-tube w/o cuff
1250-107	R1 shock-proof ® black, 1-tube, velcro cuff (adult)
1250-122	R1 shock-proof ® black, 1-tube, velcro cuff (large adult)
1250-129	R1 shock-proof ® black, 1-tube, velcro cuff (child)
1251-150	R1 shock-proof ® white, 1-tube, disinfectable one-piece cuff (adult)
1251-152	R1 shock-proof ® white, 1-tube, disinfectable one-piece cuff (large adult)
1251-154	R1 shock-proof ® white, 1-tube, disinfectable one-piece cuff (child)
1255	Set I, R1 shock-proof ® white with 3 disinfectable cuffs, 1-tube, with case latex free
1260	Set I, R1 shock-proof ® white, 3 velcro cuffs, 1-tube with storage box
1261	Set II, R1 shock-proof ® white with 3 velcro cuffs, 1-tube without box

Further items are available, please contact your Riester distributor for further information.

Please visit **riester.de** for more information