



precisa[®] N **shock-proof**

Shock-resistant aneroid sphygmomanometer

Handheld, portable aneroid sphygmomanometer with shock-resistant technology, meeting standard DIN EN ISO 81060-1. The **precisa[®] N shock-proof** aneroid sphygmomanometer is tested to be shock-resistant to a height of 1 m, providing reassurance that it will remain accurate even after impact.

- › Precision air-release valve with fine adjustment, wear-free
- › Microfilter protects air-release valve and measuring system
- › Ergonomically shaped spoon-type handle made of stainless steel
- › Linear, tilted aluminium scale Ø 63 mm, ranging from 0 to 300 mmHg (easy-to-read)
- › Latex-free
- › Metal valve base
- › Supplied in a vinyl zip case
- › Approved by the British and Irish Hypertension Society (BIHS)
- › Shock-resistant according to DIN EN ISO 81060-1
- › Shock-resistant from a height of 1 m (39.37 in)

1364-100	precisa [®] N shock-proof , 1-tube w/o cuff
1364-122	precisa [®] N shock-proof , 1-tube, velcro cuff (adult xxl)
1364-107	precisa [®] N shock-proof , 1-tube, velcro cuff (adult)
1364-129	precisa [®] N shock-proof , 1-tube, velcro cuff (child)
1364-150	precisa [®] N shock-proof , 1-tube, disinfectable one-piece cuff (adult)
1364	babyphon [®] , precisa [®] N shock-proof Set , 1-tube, 3 cuffs (child, infant and newborn)

Please visit **riester.de** for more information

