



Handheld, portable, palm-sized aneroid sphygmomanometer with centrally positioned, trigger-style, air release valve that enables left- or right-hand operation. Available in two colours – slate grey and blue – with a comprehensive selection of cuff sizes.

- Push-in tube connection at the top of the manometer for unhindered blood pressure measurement
- Latex-free bulb to inflate the cuff
- Precision pressure-release valve can be finely regulated by finger pressure without wear and tear
- Comprehensive selection of different sizes and versions of single-tube cuffs
- Integrated stethoscope version available
- Linear scale Ø 64 mm from 0 to 300 mmHg
- Maximum error tolerance ± 3 mmHg
- Microfilter for longer durability
- Supplied in a padded nylon zip bag
- Special hardened copper-beryllium membranes, virtually non-wearing
- Cuff with practical pull-through strap when purchasing a set that includes a stethoscope
- Approved by the British and Irish Hypertension Society (BIHS)

1512-100	ri-san® manometer, slate grey, w/o cuff	1517-100	ri-san® manometer, blue, w/o cuff
1512	ri-san® slate grey, velcro cuff (adult)	1517	ri-san® blue, velcro cuff (adult)
1512-122	ri-san® slate grey, velcro cuff (adult XXL)	1517-122	ri-san® blue, velcro cuff (adult XXL)
1512-129	ri-san® slate grey, velcro cuff (child)	1517-129	ri-san® blue, velcro cuff (child)
1444-141	ri-san® slate grey, 1-tube with stethoscope, velcro cuff (adult) pull through strap (D-ring)	1443-141	ri-san® blue, 1-tube with stethoscope, velcro cuff (adult) pull through strap (D-ring)

Further items are available, please contact your Riester distributor for further information.

Please visit **riester.de** for more information